






2.674
2016 Spring


Introduction to Microfluidics III

Hydrodynamics and Electrokinetics


Inside droplets: secondary flow



Chaotic Mixing in Droplets



Effect of surface roughness



Wenzel's model

- If the surface has a high free energy, roughness promotes wetting.
- If it has low free energy, roughness promotes hydrophobicity.

$$\cos \theta^* = r \cos \theta$$

$$r = \frac{\text{actual_area}}{\text{projected_area}}$$


$$\theta^* = \text{apparent_contact_angle}$$

Cassie's model

- Wettability of heterogeneous (solid+air) surfaces
- Contact angle on air fraction is 180°.


$$\cos \theta^* = -1 + \phi_s (\cos \theta + 1)$$

$$\phi_s = \text{solid_fraction_surface}$$





© sources unknown. All rights reserved. This content is excluded from our Creative Commons license. For more information, see <https://ocw.mit.edu/help/faq-fair-use>.

Cassie to Wenzel transition



A. Lafuma Nature Materials 2 457-460 (2003)

$$\cos \theta_c = \frac{\phi_s - 1}{r - \phi_s}$$



<https://en.wikipedia.org/wiki/Ultrahydrophobicity>

© sources unknown. All rights reserved. This content is excluded from our Creative Commons license. For more information, see <https://ocw.mit.edu/help/faq-fair-use>.

Electrowetting



© sources unknown. All rights reserved. This content is excluded from our Creative Commons license. For more information, see <https://ocw.mit.edu/help/faq-fair-use>.

Digital microfluidics

Electrowetting on Dielectrics (EWOD)

Droplets can be manipulated by electrowetting


Application of voltage bias between electrodes used to manipulate droplets

Complex operations can be performed

This image has been removed due to copyright restrictions.
Please see <https://mgitecetek.files.wordpress.com/2012/05/ew.jpg?w=640>.

Electrowetting

W. C. Nelson, C.-J. 'CJ' Kim / J. Adhesion Sci. Technol. 26 (2012) 1747–1771 1751




$$\frac{\text{Interfacial Force}}{\text{Length}} = C \frac{V^2}{2}$$

© Informa UK Limited. All rights reserved. This content is excluded from our Creative Commons license. For more information, see <https://ocw.mit.edu/help/faq-fair-use>.


Work of Wetting

$$W_s = \gamma_{SG} - \gamma_{SL} - \gamma_{LG} \geq 0$$




$$dW = \gamma_{LG} \cos \theta dA + \gamma_{SL} dA - \gamma_{SG} dA = 0$$

Contact Angle vs. V





$$\cos \theta = \cos \theta_0 + \frac{1}{2} \frac{\epsilon_0 \epsilon_r}{\gamma d} V^2$$

Lippmann Young Equation



Drop Manipulation


- Merge
- Split
- Move



Digital microfluidics



This image has been removed due to copyright restrictions.
Please see https://www.aaas.org/sites/default/files/content_files/1007_STM.jpg.

This image has been removed due to copyright restrictions.
Please see <https://www.technologyreview.com/s/415654/novel-chip-for-monitoring-breast-cancer/>.



© sources unknown. All rights reserved. This content is excluded from our Creative Commons license. For more information, see <https://ocw.mit.edu/help/faq-fair-use>.

Digital microfluidics



© U. Washington, UT Austin, and sources unknown. All rights reserved.
 This content is excluded from our Creative Commons license. For more
 information, see <https://ocw.mit.edu/help/faq-fair-use>.


Electrokinetics

- Electrokinetics deals with electrically driven flow of charges, particles, and fluids typically in the presence of solid-liquid interfaces


Electrokinetics

- Electrokinetics deals with electrically driven flow of charges, particles, and fluids typically in the presence of solid-liquid interfaces

Fluid Conveyor Belt



© AIP Publishing LLC. All rights reserved. This content is excluded from our Creative Commons license. For more information, see <https://ocw.mit.edu/help/faq-fair-use>.





Induced Charge Electro Osmosis (ICEO)

Electroosmotic pump,
Urbanski et al, *APL*
(2006)

Microchannel flows

- Poiseuille Flow and Couette Flow




© Cambridge University Press. All rights reserved. This content is excluded from our Creative Commons license. For more information, see <https://ocw.mit.edu/help/faq-fair-use>.

B. Kirby, *Micro- and Nanoscale Fluid Mechanics*, Cambridge Univ. Press, 2010


Poiseuille Flow

$$\nabla p = \eta \nabla^2 \vec{u}$$





Poiseuille flow in a circular tube.

$$Q = -\frac{\pi R^4}{8\eta} \cdot \frac{\partial P}{\partial z} \quad Q = \frac{\Delta p}{R_h}$$



© Cambridge University Press. All rights reserved. This content is excluded from our Creative Commons license. For more information, see <https://ocw.mit.edu/help/faq-fair-use>.

Hydrodynamic Force Assisted Sorting



© sources unknown. All rights reserved. This content is excluded from our Creative Commons license. For more information, see <https://ocw.mit.edu/help/faq-fair-use>.

Couette Flow


Navier-Stokes

$$\rho \frac{\partial \vec{u}}{\partial t} + \rho \vec{u} \cdot \nabla \vec{u} = -\nabla p + \eta \nabla^2 \vec{u} \quad \longrightarrow \quad 0 = \eta \frac{\partial^2 u}{\partial y^2}$$


This image has been removed due to copyright restrictions.
Please see <http://www.kirbyresearch.com/images/etc/textbook/mae28.jpg>.

Electrophoresis and Electroosmosis

- Electrophoresis: Movement of charged particles under the influence of electric field



- Electroosmosis: Flow of fluid under influence of electric field




Origins of surface charge

Chemical origins:

- Dissociation of surface groups
- Adsorption of charged species

Physical origins:



- Induced charge



This image has been removed due to copyright restrictions. Please see https://commons.wikimedia.org/wiki/File:Double_Layer.png.

© sources unknown. All rights reserved. This content is excluded from our Creative Commons license. For more information, see <https://ocw.mit.edu/help/faq-fair-use>.

Electrical double layer



B. Kirby, Micro- and Nanoscale Fluid Mechanics, Cambridge Univ. Press, 2010

© Cambridge University Press. All rights reserved. This content is excluded from our Creative Commons license. For more information, see <https://ocw.mit.edu/help/faq-fair-use>.

Electroosmosis

Assumptions:

$$u = u(y)$$

$$v = w = 0$$

$$E_{xy} = \text{constant}$$

$$\rho = \text{constant}$$

This image has been removed due to copyright restrictions. Please see <https://www.kirbyresearch.com/images/etc/textbook/mae28.jpg>.

- Macroscale Navier-Stokes eqn.

- N-S equation with negligible pressure drop, steady flow $u(y)$, and surface charge

- Poisson eqn. (electrostatics)

B. Kirby, Micro- and Nanoscale Fluid Mechanics for Engineers: Transport in Microfluidic Devices <http://www.kirbyresearch.com/index.cfm/wrap/textbook/microfluidicsnanofluidics.html>.

Electrical double layer

Gouy-Chapman Model

- Boltzmann distribution


$$n_{\pm} = n_{bulk} \exp\left(\frac{\mp e\phi}{kT}\right)$$

- Poisson-Boltzmann equation

$$\nabla^2 \phi = \frac{-\rho_E}{\epsilon} = -\frac{e(n_+ - n_-)}{\epsilon}$$


- Debye length

$$\lambda_D = \sqrt{\frac{\epsilon kT}{2n_{bulk} e^2}} \quad (1-100 \text{ nm})$$




Electrical double layer

- Debye length
 - 1 KCl \rightarrow \sim 0.3 nmh
 - 100 mM KClh \rightarrow 1 nmh
 - 1 m KClh \rightarrow 10 nmh
 - Water \rightarrow 1 μ m



- Zeta potential
 - Typically 0-100 mVh
 - Determined by surface charge h and ionic concentration



Electroosmosish

This image has been removed due to copyright restrictions. Please see <https://www.kirbyresearch.com/images/etc/textbook/mae28.jpg>.

Assumptions:h

$$u = u(y)$$

$$v = w = 0$$

$$E_{xy} = \text{constant}$$

$$\rho = \text{constant}$$

$$\eta \frac{d^2 u}{dy^2} = \varepsilon \frac{d^2 \phi}{dy^2} E_x$$

$$\eta u = \varepsilon E_x \phi + C_1 y + C_2$$

BCs: h

Velocity is finite as $y \rightarrow \infty$


$$u(0) = 0$$

$$u = \frac{\varepsilon E_x}{\eta} (\phi - \phi_0)$$

$$u_{eo} = -\frac{\varepsilon \zeta}{\eta} E_x$$

(Helmholtz-Smoluchowski equation)h


Vorticity Generation



Vorticity generation and cancellation in EDLs.

B. Kirby, Micro- and Nanoscale Fluid Mechanics for Engineers: Transport in Microfluidic Devices
<http://www.kirbyresearch.com/index.cfm/wrap/textbook/microfluidicsnanofluidics.html>


Zeta potentials of glass and PDMS




Kirby & Hasselbrink, *Electrophoresis* (2004)

© John Wiley & Sons, Inc. All rights reserved. This content is excluded from our Creative Commons license. For more information, see <https://ocw.mit.edu/help/faq-fair-use>.

Electroosmotic pump




Electroosmotic pump



→ Electroosmotic pumping is more effective at small length scales!

Electrophoretic mobility




- Presence of counterions affects mobility

Electrophoretic mobility

- Mobility (μ): Particle velocity under unit electric field

$$\vec{V}_{ep} = \mu \vec{E}$$



$$F_{drag} = 6\pi\eta r V$$


$$\vec{V}_{ep} = \frac{q\vec{E}}{6\pi\eta r}$$

$$\mu_{EF} = \left(\frac{q}{6\pi\eta r} \right)$$

$$F_{electric} = qE$$

(Hückel equation)


Capillary Electrophoresis (CE)



This image has been removed due to copyright restrictions.
Please see the image on Page 10 of
<http://nanoparticles.org/pdf/noh.pdf>.

Electrophoretic/electroosmotic mobilities


- Small monovalent ions (Na^+ , K^+ , Cl^- , etc.): $\sim 10^{-7} \text{ m}^2/\text{v.s}$
- Cells: $\sim 10^{-8} \text{ m}^2/\text{v.s}$
- In PDMS, glass: $\sim 10^{-7}$ to $10^{-8} \text{ m}^2/\text{V.s}$
- DNA: $4 \times 10^{-8} \text{ m}^2/\text{V.s}$ or lower (in gels)



DNA mobility in TAE buffer
Salieb-Beugelaar et al., *Lab Chip* 9, 2508 (2009)

Image courtesy of U.S. National Library of Medicine.

Electrophoretic separation




Separation based on differential migration velocities

Migration velocities depend on

- Minimize electroosmotic flow, PVP (polyvinylpyrrolidone)
- Size and charge of particles
- Interaction with gels/pores, if any

$$\mu_{EF} = \left(\frac{q}{6\pi\eta r} \right)$$

Electrophoretic separation



$$\Delta_{1-2} = (V_2 - V_1)t = (\mu_{ep,2} - \mu_{ep,1})Et$$


$$\Delta_1 \sim \sqrt{D_1 t}$$

$$\frac{\Delta_{1-2}}{\Delta_1} = \frac{(\mu_{ep,2} - \mu_{ep,1})Et}{\sqrt{D_1 t}} \propto E\sqrt{t}$$


Electrophoresis in microfluidics

- Advantages:
 - Higher electric fields easily possible
 - Shorter time-to-result
 - Integration with other processing steps
 - Low sample requirements

Step 1: Injection of buffer solution




Step 2: Injection of sample solutions




39

Step 3: Start of electrophoresis



40

Step 4: EP Separation



41

Capillary Electrophoresis (CE)

This image has been removed due to copyright restrictions.
Please see the image on Page 9 of <http://nanoparticles.org/pdf/noh.pdf>.

This week

- Examine the motion of microparticles under electric fields
- Study the effect of surface treatment on electrokinetic flow
- Separate and analyze molecules using gel electrophoresis

MIT OpenCourseWare
<https://ocw.mit.edu>

2.674 / 2.675 Micro/Nano Engineering Laboratory
Spring 2016

For information about citing these materials or our Terms of Use, visit: <https://ocw.mit.edu/terms>.