

**6.542J, 24.966J, HST.712J LABORATORY ON THE PHYSIOLOGY,
ACOUSTICS, AND PERCEPTION OF SPEECH
Fall 2005**

Kenneth N. Stevens
Stefanie Shattuck-Hufnagel

Joseph S. Perkell

Lab 1 Handout

09/08/05

Using Klatttools on UNIX system

- 1. To log in**
Type *labc*
Type *password*
Type *startx*
Click on *xterm*
Cd *group1, group2, group3, ...* (depending on your group)
- 2. To make a spectrogram of a file xxx.wav**
Type *lspecto xxx*
This creates a ps file *xxx.ps*
To print the spectrogram type *lpr xxx.ps*
- 3. To examine waveforms and spectra**
Type *xkl xxx*
This will display a waveform (at top), an expanded waveform (bottom left), and a spectrum at the location of the cursor (middle left).
Parameters of spectrum can be changed by clicking on *spectrum*.

To make a spectrogram

Type *i*

To save a spectrum or waveform

Click on file and follow instructions.

(a) open a ps file *yyy.ps* (b) save to ps file (c) close ps file

To print, type *lpr yyy.ps*

- 4. To edit a waveform**

Set start of selected region using cursor, type *w*.

Set end of selected region using cursor, type *e*.

Save *w - e* to *.wav* file by typing *o*.